


TransWin

For:

- REA ScanCheck
- REA Check 3

The TransWin program is intended to transfer and administrate verification protocols of REA Bar Code Verifiers



Functions:

- Selectable language English/German
- Display of verification protocol and bar code reflection profile on the monitor
- Facility to backup the verification protocols and the bar code reflection profile
- Powerful search function for open and saved verification protocols by using lists of terms and logical operators
- Easy to use print function that is able to add company information and comments.
- Option to save verification protocols as text which can be used with a word processor

System requirements

- Windows 3.1 / 3.11 / 95/NT
- One free serial port
- REA ScanCheck or REA Check 3 verifier.

The verification protocols are transferred by using the print function via a standard serial port to a Windows PC.

Each verification protocol is shown in a separate window. The reflection profile of the REA ScanCheck has also a separate window.

For backup purposes the verification protocols with reflection profile can be saved to a file.

TransWin is able to display directly the received or saved verification protocols.

To identify the user it is possible to automatically add a company description to each verification protocol.

A search function offers the option to look through the saved or in memory protocols. It is possible to search for predefined terms or for directly typed in terms. The predefined terms can be combined by using the logical functions AND or OR.

REA Elektronik GmbH

Rationalisierung · Entwicklung · Automation
Elektronische Geräte und Systeme
Teichwiesenstraße 1
D-64367 Mühlthal-Waschenbach
Tel.: + 49 6154 638-0
Fax: + 49 6154 638-195
email: reainfo@rea.de
Internet: http://www.rea.de